

Excellent Illumination Power In A Space-Saving, Affordable Unit



Introducing The CLH-SC Halogen Light Source – Compact, Powerful, And Reasonably Priced

Extraordinarily compact yet offering light intensity powerful enough for monitor observation, the CLH-SC is the ideal light source for the OTV-SC Video System. Space efficiency, and high power make the CLH-SC an excellent choice for monitor observation on an outpatient basis.



Powerful Illumination

In spite of its diminutive size, the CLH-SC's 150 W halogen lamp generates more than enough light intensity to fully illuminate an observation site and its surroundings. As a result, it can be clearly reproduced on your monitor.

Space-Efficient, Compact Design

With a width of only 220 mm and a depth of 290 mm, the CLH-SC can easily be installed in confined spaces. Since it has the same dimensions as the OTV-SC Video System, it can be placed beneath or on a top the OTV-SC to achieve maximum space efficiency in your office. Small size combined with a lightweight design also makes the CLH-SC easy to move whenever necessary.

Item	Specification	
Power Supply	Voltage	120V AC 220-240V AC
	Voltage fluctuation	Within $\pm 10\%$
	Frequency	50/60 Hz
Size	Dimensions	220(W) x 74(H) x 290(D) mm
	Weight	5 kg
Illumination	Lamp	Halogen lamp 150 W
	Lamp Life	Approx. 50 hours on continuous use
	Brightness control	Light-path diaphragm by using the dial on the front panel
	Lamp cooling method	Forced-air cooling
Classification	Type of protection against electric shock	Class I
	Degree of protection against electric shock	TYPE BF
	Degree of protection against explosion	Product should be kept away from flammable gases
Operating Environment	Ambient temperature	10 to 40° C (50 to 104° F)
	Relative humidity	30 to 85%
	Air pressure	700 to 1,060 hPa



Specification, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS[®]

OLYMPUS OPTICAL CO., LTD.

2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
OLYMPUS WINTER & IBE GMBH
 Kuehnstrasse 61, 22045 Hamburg, Germany

OLYMPUS OPTICAL CO. (EUROPA) GMBH.

Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany

OLYMPUS AMERICA INC.

2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.

OLYMPUS LATIN AMERICA, INC.

6100 Blue Lagoon Drive, Suite 390 Miami, Florida 33126-2087, U.S.A

OLYMPUS KEYMED

Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.

491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373

OLYMPUS HONG KONG AND CHINA LIMITED.

Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong

OLYMPUS BEIJING REPRESENTATIVE OFFICE

Rm.818, South Tower, Beijing Kerry Centre, No.1 Guanghua Road, Chaoyang District, Beijing 100020, China

OLYMPUS MOSCOW LIMITED LIABILITY COMPANY

117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia

OLYMPUS AUSTRALIA PTY. LTD.

31 Gilby Road, Mount Waverley, VIC., 3149, Australia

OLYMPUS BUSINESS AREAS

IMAGING

MEDICAL

INDUSTRIAL

www.olympus.com